
HANDHABUNGSHINWEIS  
INSTRUCTIONS FOR USE  
INSTRUCTIONS DE SERVICE  
ISTRUZIONI PER L'USO  
INSTRUCCIONES DE USO

Monopolare MIC Elektroden  
Monopolar MIS electrodes  
Electrodes monopolaires CMI  
Elettrodi monopolari CMI  
Electrodos monopolares CMI

85900000 - 85902403

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**Monopolar electrodes for MIS**

**Attention**

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

**Intended Use**

The monopolar electrodes are intended for use in minimally invasive and especially laparoscopic surgical procedures. The device is to be passed through a 5.5 mm laparoscopic cannula.  
The monopolar electrodes micro needle electrodes are intended to be used for monopolar coagulation and cutting respectively. The monopolar hysterectomy loop electrode is intended to be used for removal of the uterus in supracervical hysterectomies.

The instrument has to be connected by a suitable monopolar cable to the monopolar output of an electrosurgical generator.

As indicated monopolar coagulation- or cutting current may be selectively applied to the tissue.  
The maximum output voltage of the generator must not exceed 2000 Vp.

**Suitable Cables:**

Bissinger Monopolar Cable REF 801 00xxx.

**Attention:** Instruments for electrocautery should be used only by persons who have been specially trained in the use of such instruments.

**Contraindications**

Incidents which have been reported in connection with the use of electrosurgical systems:

Unintended activation with resulting tissue injury on the wrong spot and/ or damage to the equipment.  
Fire in connection with surgical drapes and other inflammable materials.  
Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.  
Explosions caused by sparks in the proximity of inflammable gases.  
Perforation of organs. Sudden severe bleedings.  
Bipolar forceps have proved inefficient for tubular sterilisation or coagulation in the context of sterilisation and should therefore not be used for this purpose.

**Use and safety instructions**

The non-observance of the present use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.

- Before initial use and any other use, all instruments have to be completely cleaned, disinfected and sterilised and their functioning has to be checked.
  - It is very important to check each surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas, such as blades, tips, notches, locking and blocking devices as well as all mobile parts, insulations and ceramic elements have to be checked carefully.
  - Never use damaged instruments.
  - Never use the instruments in the presence of flammable or explosive substances.
  - The instrument may not be laid down on the patient.
  - Coagulation should only be performed if the contact surfaces are visible and ensure a good contact to the tissue selected for coagulation. Do not touch any other metallic instruments, trocar sleeves, optics or things like that during use.
  - Do not use in patients with pacers or other active implants without seeking the corresponding technical advice first.
  - Ensure correct application of the neutral electrode on the patient; otherwise, there is a danger of burns.
- Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.

**Assembly**

Follow the pictogram PIC-859 for assembly and disassembly of the instrument.

**Special instructions for the hysterectomy loop electrode**

- Insert the instrument through the cannula with the loop pulled inside the shaft.
- Thoroughly check correct placement of the loop around the uterus prior to activation of the electrical current; the entire tissue surrounded by the loop will be cut.
- After correct placement, maintain tension of the loop to avoid uncontrolled slip of the loop. Maintain some safety distance to the ligature.
- During electrosurgical cutting, pull the loop continuously inside the shaft. Cutting will start at the uninsulated part of the loop wire.
- The loop wire is a replaceable, sterilisable product for single use that should be replaced following each procedure.

**Reprocessing**

Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling.  
Instruments for electrocautery are by nature subject to increased wear depending on the type and time of use.

**Preparation and transport**

Remove coarse dirt from the instruments immediately after each use. Do not use fixation agents or hot water (>40°C). Storage and transport of the instruments to the reprocessing location must be ensured in a sealed container.

**Machine reprocessing**

**Cleaning**

Place the instruments in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.

1. Pre-rinse for 1 min. with cold water
2. Discharging
3. Pre-rinse for 3 min. with cold water
4. Discharging
5. Wash for 5 min. at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent.
6. Discharging
7. Neutralise for 3 min. with warm tap water (>40°C) and a neutralising agent.
8. Discharging
9. Rinse for 2 min. with warm tap water (>40°C).
10. Discharging

**Disinfection**

Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

**Drying**

Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.  
If necessary, manual drying may additionally be carried out using a lintfree cloth. Dry cavities by blowing with sterile compressed air.

**Manual reprocessing**

**Ultrasonic Pre-Cleaning**

1. The instrument must be inserted in an ultrasonic bath with 0.5% enzymatic cleaning detergent. Ultrasonic must be applied for 15 minutes at 40°C/104°F.
2. Remove the instrument and rinse completely with cold water to remove the cleaning detergent.

**Cleaning**

Prepare a cleaning bath according to the manufacturer's instructions.

1. Never use products with cold tap water (<40°C) until all visible accumulations of dirt have been removed. Remove stuck dirt by using a soft brush.
2. Place products in the prepared cleaning bath so that they are completely submerged. Observe residence time according to the manufacturer's instructions.
3. Clean the instrument in the bath manually using a soft brush. Brush all surfaces several times.
4. Rinse the products thoroughly with DI water to remove the cleaning agents without residue.

**Disinfection**

Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer. Place the instruments in the disinfectant bath and observe the specified residence time. Rinse the products thoroughly with fully demineralised water to remove the disinfectant without residue.

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SMP GmbH # 01707011901 (limpieza mecánica)
MDS GmbH # 135196-10 (limpieza/desinfección manual)
MDS GmbH Testbericht 084183-10 (esterilización)