HANDHABUNGSHINWEIS INSTRUCTIONS FOR USE INSTRUCTIONS DE SERVICE ISTRUZIONI PER L'USO INSTRUCCIONES DE USO

Non-Stick MONOPOLARES SAUG-/KOAG.INSTRUMENT

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Reinigungsmittel (maschinell):

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Reinigungsmittel (manuell):

Enzol Enzym. Detergent, Johnson&Johnson

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ENGLISH Monopolar Suction-Coagulator

Attention

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

Intended Use

The monopolar suction-coagulator is intended to be used for the aspiration of blood, liquids and smoke and simultaneous or subsequent monopolar coagulation of tissue during surgical procedures.
The instrument has to be connected by a suitable monopolar cable to the monopolar output of an electrosurgical generator and must only be used with monopolar coagulation current.
The maximum output voltage of the generator must not exceed 3000 V_p.

Suitable Cables:

Bissinger Monopolar Cable REF 801 00xxx.

Caution: Monopolar instruments should be used only by individuals who are trained and licensed to use such devices.

Examples of such training and experience include: Training through qualified residency program, surgical skills workshops, training programs offered by equipment manufacturers or preceptorship/surgical assistant training.

Use and safety instructions

- All instruments have to be completely cleaned, disinfected and sterilised before initial use and any subsequent use.
- It is very important to check each surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use.
- Never use damaged instruments.
- Never use on patients with pacemakers or other active implants without prior consultation of a specialist.
- Never use the instruments in the presence of flammable or explosive substances.
- Ensure correct application of the neutral electrode on the patient. Otherwise there is a risk of burns.
- The instrument may not be laid down on the patient.
- Activation should only be performed if the contact surfaces are visible. Do not touch any other metallic instruments during coagulation.
- Activate in short sequences, avoid uninterrupted activation >15 sec.
- Do not use in a Spray-/Fulguration mode.
- Use non-conductive liquid when used together with irrigation instruments.

Reprocessing

Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling.
Instruments for electrosurgery are by nature subject to increased wear depending on the type and time of use.

Preparation and transport

Remove coarse dirt from the instruments immediately after each use (within 2 h). Do not use fixation agents or hot water (>40°C) as this may result in the fixation of residues and could reduce the cleaning success.
Storage and transport of the instruments to the reprocessing location must be ensured in a sealed container to avoid any damage to the instruments and any contamination of the environment.

Manual Pre-Cleaning

1. Immerse the instrument in cold water for 5 minutes.
2. Brush the instruments under cold water until all visible impurities are removed.
3. Suction tube, suction control holes and safety valves must additionally be cleaned with the Bissinger cleaning brush (REF 89800000) at least six times.

4. All lumen must be flushed with the water pistol for 10 seconds or using a syringe with at least 50 ml water several times, thereby holding suction control holes closed. Repeat procedure until no impurities are visible in the flushing solution.

Suitable cleaning brush:

Bissinger Cleaning Brush REF 89800000

Machine reprocessing

Cleaning

Place the instruments in a basket on the insert module of the washing machine and start the cleaning process.

1. Pre-rinse for 4 min. with cold water
2. Discharging
3. Wash for 5 min. at 50°C/122°F with a 0.2% alkaline cleaning agent (pH > 10).
4. Discharging
5. Neutralise for 3 min. with warm tap water (>40°C/104°F) and a neutralising agent.
6. Discharging
7. Rinse 2 min. with cold tap water.
8. Discharging

Disinfection

Machine operated thermal disinfection at 90°C has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

Drying

Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.
If necessary, manual drying may additionally be carried out using a lint free cloth. Dry all cavities of the instruments by blowing with sterile compressed air.

Manual reprocessing

Ultrasonic Pre-Cleaning

1. The instrument must be inserted in an ultrasonic bath with 0.5% enzymatic cleaning detergent. Ultrasound must be applied for 15 minutes at 40°C/104°F.
2. Remove the instrument and rinse completely with cold water to remove the cleaning detergent.

Cleaning

1. Prepare a cleaning bath according to the manufacturer's instructions.
2. Rinse products with cold tap water (<40°C/104°F) until all visible accumulations of dirt have been removed. Remove stuck dirt by using a soft brush. All cavities, lumen and holes must additionally be cleaned with the Bissinger cleaning brush (REF 89800000) at least six times.
3. Place products in the prepared cleaning bath so that they are completely submerged. Observe residence time according to the manufacturer's instructions.
4. Clean the instrument in the bath manually using a soft brush. All surfaces have to be brushed several times.
5. **The following steps only apply to channels and the insides of tubes:** The Bissinger cleaning brush has to be pushed in and out of the suction tube, the suction control holes and the safety valves at least six times. Rinse the tubes with distilled water and repeat the procedure, thereby holding suction control holes closed.
6. Rinse the products thoroughly with running tap water to remove the cleaning agents without residue.

Disinfection

Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer.
Place the instruments in the disinfectant bath and observe the specified residence time.
Rinse the products thoroughly with fully demineralised water to remove the disinfectant without residue.

Drying

Manual drying is carried out using a lint free cloth and, in particular, for drying cavities and channels, sterile compressed air.

Functional test and packaging

Perform visual inspection for cleanliness; if required, perform an assembly and functional test according to the operating instructions.
If necessary, repeat the reprocessing process until the instrument is optically clean.
Packaging has to comply with international standards for packaging for sterilised instruments.

The noble metal tips of the Non-Stick electrodes may tarnish like silver. This has no impact on the functionality of the instrument. The original gloss can be re-established with a silver cleaning cloth.

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- Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) in consideration of the respective national requirements.
- 3 pre-vacuum phases with a pressure of at least 60 mbar.
 - Heating up to a sterilisation temperature of min. 132°C and max. 137°C
 - Shortest exposure time: 3 min.
 - Drying time: at least 10 min.

Storage

Sterilised instruments have to be stored in a dry, clean and dust-free area at moderate temperatures from 5°C to 40°C.

Repairs

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. For any question on these matters, please contact either the manufacturer or your medico-technical department.

Attention: Defect products must pass the complete reprocessing process before being returned for repair.

Information on the validation of the reconditioning

The following testing instructions, materials and equipment have been used for validation:

Cleaning agents (for machine use):

- Neodisher FA by Dr. Weigert (alkaline)
- Endozyme by Ruhof (enzymatic)

Cleaning agents (manual cleaning):

- Enzol Enzym, detergent by Johnson&Johnson

Disinfectants (manual disinfection):

- Cidex OPA, Johnson&Johnson

Neutrising agent:

- Neodisher Z by Dr. Weigert
- **Cleaning and disinfection device:**
- Miele G 7736 CD
- Miele insert module E 327-06
- Miele MIS module E 450

For details, see report.

- SMP GmbH # 01707011901 (machine cleaning)
- MDS GmbH # 135196-10 (man. cleaning/disinfection)
- Nelson Labs # 200432706-02 (sterilisation)
- MDS GmbH Testbericht 084183-10 (sterilisation)

If the chemicals and machines described above are not available, the user has to validate the use process accordingly.

Handling

During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care.

This applies particularly to blades, fine tips and other sensitive areas.

Warranty

Günter Bissinger Medizintechnik GmbH exclusively supplies tested and faultless products to the customers. All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which, compared to the original product, have been modified, misused or handled or used in an inexpert way.

Explanation of symbols

	LOT-Number
	Caution: Non-sterile product
	Reference number
	Attention !
	Refer to instructions for use

CE-Mark and registration number of the Notified Body DQS Medizinprodukte GmbH August-Schanz-Straße 21 60433 Frankfurt, Germany

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SMP GmbH # 01707011901 (limpieza mecánica)
MDS GmbH # 135196-10 (limpieza/desinfección manual)
MDS GmbH Testbericht 084183-10 (esterilización)