## MICRODOME

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Bipolar coagulation electrode  
for minimal access procedures

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
## Microdome

### Bipolar coagulation electrode for

#### minimal access procedures

Surgery is more and more frequently effected through minimal access points, either through trocar tubes or through natural body orifices. This requires instruments which, on the one hand, have the necessary stability and, on the other hand, are sufficiently fine to ensure that the surgical areas are always visible despite of the small access.

In cases where bipolar forceps cannot be used e.g. for bipolar coagulation of diffusely bleeding areas in the nose, a new generation of electrodes have been developed with a shaft diameter of only two mm. Microdome electrodes enable the surgeon to cauterise bleeding in an effective and safe way even in the smallest cavities.

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
#### Microdome Bipolar coagulation electrode for minimal access procedures

	Bajonett Bayonet	gewinkelt 20° angled 20°	biegbar malleable
13 cm	Art. No. 856 20410	Art. No. 856 20420	—
14 cm	—	—	Art. No. 856 20430

Bild/ figure 2

The two electrical poles form a divided cupola on the working tip, (figure 1) so that when the bipolar coagulation current is activated, a coagulation area with a diameter of approximately 2.5 mm and a depth of approximately 0.5 mm is clearly defined (figure 2). The desired coagulation result is obtained at a very low generator power of about 10 watts so thermal damage to surrounding tissue is kept to a minimum.

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The rounded cupola tip aids insertion into the surgical area and by using bipolar technology, the application of a neutral electrode is not required.

The ergonomic handle is equipped with a standard flat plug connection that can be used with common bipolar cables to nearly all high-frequency generators. The instrument can be cleaned easily as it does not have to be dismantled, and is completely reusable; as a result the costs per use are low.

The Microdome coagulation electrodes are available in angular or bayonet forms. The design maintains the visibility of the surgical area by keeping the hand holding the instrument out of the line of sight and so safe working is ensured. Furthermore, a malleable electrode model is available which surgeons can bend themselves, so that individual adaptation of the instrument to particular anatomical features is possible (figure 3).

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**Bissinger – Focus on better solutions**

Bipolar and monopolar instruments  
 Accessories for electrosurgery  
 Surgical instruments  
 Development of instruments  
 to customer specifications  
 Service and repair

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**Medizintechnik GmbH**