


89505210

CeraNeedle Micro Dissections Needle Electrode



Technical specifications

Dimensions and Weights

Total length	60 mm
Dimensions	mm
Diameter	0,5 mm

Basic Specifications

Shaft diameter	2,4 mm
Electrode shape	Needle

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 4,3 kVp
Packaging unit	1
GTIN 14	04250418729215