

89505212

## CeraNeedle Micro Dissections Needle Electrode



### Technical specifications

#### Dimensions and Weights

Total length	70 mm
Dimensions	mm
Diameter	0,5 mm

#### Basic Specifications

Shaft diameter	2,4 mm
Electrode shape	Needle

#### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

#### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

#### Additional Product Information

Accessories rated voltage	max. 4,3 kVp
Packaging unit	1
GTIN 14	04250418729222