


**89725110**

Monopolar Forceps, straight



## Technical specifications

### Dimensions and Weights

Total length	250 mm
Total length in inches	9,84 "
Tip width	2,1 mm

### Basic Specifications

Forceps shape	straight
---------------	----------

### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

### Additional Product Information

Accessories rated voltage	max. 2 kVp
Packaging unit	1
GTIN 14	04250418730198